


What Can EtCO₂ Do For You?

It's as easy as **PQRST** – remember this waveform to help you use EtCO₂ effectively


R OSC

Return of spontaneous circulation is readily evident following a successful defibrillation shock; EtCO₂ rises rapidly to over 40 mmHg.




QUALITY

CPR quality is reflected in the EtCO₂ waveform. If it falls off, time to switch rescuers! EtCO₂ will also alert you to hyperventilation.



STRATEGY


EtCO₂ may indicate some Hs and Ts. For example, when hypoventilation occurs, a gradual increase in EtCO₂ is displayed on the capnogram.



R

PLACEMENT

Be sure the endotracheal tube is in the right place and doesn't get dislodged.



TERMINATION

- An end-tidal carbon dioxide level of 10 mmHg or less, measured 20 minutes after the initiation of advanced cardiac life support, accurately predicts death in patients with cardiac arrest associated with electrical activity but no pulse.
 - Arrest survivors' EtCO₂ averaged 32.8 ±7.4 mmHg 20 minutes into the arrest.
 - Non-survivors' EtCO₂ averaged 4.4 ±2.9 mmHg at 20 minutes.
- Source: Levine, et. al. *New Engl J Med*, 337:301-6, 1997

P

Q

S

T



The X Series® monitor/defibrillator with EtCO₂ monitoring provides valuable clinical information at a glance.

REVERSIBLE CAUSES OF CARDIAC ARREST

- Hypovolemia
- Hypoxia
- Hydrogen ion (acidosis)
- Hyperkalemia
- Hypokalemia
- Hypothermia
- Hypoglycemia
- Hyperglycemia
- Tablets or Toxins
- Tamponade, cardiac
- Tension pneumothorax
- Thromboembolism
- Thrombosis

Learn more at www.zoll.com or call 800-804-4356.